



**Raven Labs**  
*a division of Mesa Labs*  
8607 Park Drive  
1-800-728-5702 or 1-402-593-0781  
[info@ravenlabs.com](mailto:info@ravenlabs.com)

## Instructions for Use


### Autoclave Tape

1. Affix tape to item to be sterilized. Tape may be used to secure packages or solely as an external process indicator.
2. Sterilize items per in-house procedures.
3. Verify tape has changed color. Tape is intended to identify items having been exposed to a sterilization process.

Interpretation:



**Raven Labs**  
*a division of Mesa Laboratories, Inc.*  
8607 Park Drive, Omaha, NE 68127 USA  
[www.ravenlabs.com](http://www.ravenlabs.com)

   
EC Authorized Representative  
WMDE  
Kruisstraat 108  
6461 HC Kerkrade  
The Netherlands  
Tel: +31-(0)45-535.0625  
**0344**