



Bowie-Dick Type Tests

Raven Labs offers three types of air removal Bowie-Dick type tests intended to evaluate the performance of the air removal (vacuum) system of pre-vacuum equipped steam sterilizers. All of our Bowie-Dick tests are **lead-free** and offer an easy to interpret, distinct, non-reversible color change.

ProChem B-D

Pre-assembled, disposable test pack



Catalog # CI-BDTLF-20

Box of 20 Tests

Catalog # CI-BDTLF

Case of 4 boxes (80 Tests)

3 year expiry

- Distinct, non-reversible pink to black color change
- Meets AAMI, EN 867-3 and ISO 11140 air removal

UNEXPOSED



FAIL



PASS

After exposure, open pack and remove the indicator sheet for your records.

ProChem B-D (EN Compliant)

Pre-assembled, disposable test pack

International Use Only

Catalog # CI-BD867-4-20

Box of 20 Tests

Catalog # CI-BD867-4

Case of 4 boxes (80 Tests)

3 year expiry

- Distinct, non-reversible yellow to black color change
- Complies with EN Standard 867-4

After exposure, open pack and remove the indicator sheet for your records.



PASS

UNEXPOSED

FAIL

ProChem BD I 10

Individual Test Sheets

Industrial/International Use Only

- Distinct, non-reversible blue to black color change
- Meets AAMI, EN 867-3 and ISO 11140 air removal

- Meets AAMI, EN 867-3 and ISO 11140 air removal

UNEXPOSED

FAIL



PASS

Catalog # CI-BDI 10

100 Sheets/package

3 year expiry

Place sheet within stack of surgical towels that are assembled according to AAMI guidelines.

M E S A

L A B S